

Precision meets efficiency









State of the art aneroid palm style sphygmomanometer available as 1- and 2-tube version with disinfectable latex-free one-piece cuffs.

- > Fast inflation of the cuff due to new ergonomic spoon design
- Disinfectable and washable one-piece cuffs in different sizes
- > Comfort air release valve, wearfree, with fine adjustment.
- Large manometer scale (ø 63 mm) and optimized scale and pointer design for better readability of values.
- > High-polished manometer surface, easy to disinfect.
- Maximum error tolerance ±3 mmHg.

- > Microfilter protects air release valve and measuring system.
- Available as 1-tube and 2-tube model with plug-in connection.
- > Completely latex-free version with shock-protecting polycarbonate
- > Delivery includes one-piece cuff and vinyl zipper case.
- Approved by the BIHS validation

1370-150	e-mega® white, 1-tube, disinfectable one-piece cuff (adult)
1375-150	e-mega® black, 1-tube, disinfectable one-piece cuff (adult)
1371-151	e-mega® white, 2-tube, disinfectable one-piece cuff (adult)
1376-151	e-mega® black, 2-tube, disinfectable one-piece cuff (adult)

^{*} Further items are available, please contact your Riester distributor for further information.

Please visit riester.de for more information



79640 Rev. A 2023-10-19